


## MASTconnect Software Modules

The MASTconnect software controls the system via user inputs. The sample pilots must be each identified and enabled. The user is free to choose how often, what type, and where samples are taken and directed. The software provides error notifications and records the final data results in real time as sample tests are finished executing. The MAST software layout is user friendly and helps you achieve the variety of samples needed to maximize a bioreactors output.



### FEATURES

- **Simple layout and directives makes software simple to configure**
- **Provides information on system operation and status routinely**
- **Collects and ensures data integrity by all analytical devices**
- **Modular design allows easy reconfiguration of sample pilots and destination types.**
- **Saves and executes sample schedules**
- **Prompts user actions that will maintain system functionality**

### SPECIFICATIONS

#### Measurements:

- **NA**

#### Utilities and/or Components Required:

- **Connection to MAST System**