

Single-Use Bioreactor – Probe Sheath Adapter

The SP200 Probe Sheath Adapter is designed to connect an SP200 to a single-use bioreactor using the common aseptically connected probe sheath that is used to introduce a typical pH or DO probe.



FEATURES

- Allows an aseptic connection between the SP200 and a single use bioreactor
- Utilizes the same interface as a common pH or DO probe

SPECIFICATIONS

Product Characterization	Specification
Width	2.4 in (6.1 cm)
Height	5.5 in (14.0 cm)
Depth	11.3 in (28.7 cm)
Weight	1.1 lbs (0.5 kg)
Required Components	Sample Pilot SP200 Aseptic Probe Sheath