

## Sanitant Supply Station

This 20 L drum contains the sanitant that is flushed through the system in between each sample drawing event. By using 70% Isopropyl Alcohol or 70% Ethanol, the system is thoroughly flushed and sanitized each time sample is sent through the lines. The sanitant station is equipped with level detection alarms, which will notify the user on the software that the sanitant must be replaced.



### FEATURES

- Automatically cleans the system throughout operation
- Pumps sanitant at high enough pressure to clear out possible blockage and residue
- Utilizes level detection alarms to notify when sanitant is running low
- Simple to setup and refill

### SPECIFICATIONS

Product Characterization	Specification
Max Width	17.2 in (43.7 cm)
Max Height	21.4 in (54.4 cm)
Max Depth	12.5 in (31.8 cm)
Weight (empty)	8.0 lbs (3.6 kg)
Required Components	Connections to MAST System
	70 % IPA or 70% EtOH